



Decoding the Lens


A Clinical Guide to Lensometry, Refraction, and Optical Measurement



The Physics of Correcting Refractive Errors




Minus Lenses: Diverge light to compensate for a longer eye.



Plus Lenses: Converge light to compensate for a shorter eye (or flat cornea).


Categorizing Primary Lens Properties

Sphere



Corrects Myopia or Hyperopia.
Carries equal refractive power in
all meridians.

Toric




Corrects Astigmatism (where light
focuses on the retina in one plane,
but not another). Contains
multiple focal planes and requires
exact axis alignment.

Multifocal



Corrects Presbyopia (the natural
loss of lens accommodation).
Houses multiple focal points for
distance and near vision.


Instruments of Precision



AUTOMATED LENSOMETER

Advantage: High efficiency, prints exact Rx, ideal for high-volume clinic flow.

Limitation: Can struggle with heavily scratched lenses or complex prism measurements.



MANUAL LENSOMETER

Advantage: The clinical gold standard. Required mastery for COT and COMT JCAHPO certifications.

Limitation: Essential for troubleshooting, verifying ground-in prism, calculating slab-off, and precise progressive mark reading.

Anatomy of the Manual Lensometer

1 The Eyepiece:

The Eyepiece is used by the microscope to spectacles, into the height, the wheel and on the inversion.

2 The Axis Wheel:


The Axis wheel is a small, knurled wheel and the see observable. The spectacles and graph its on an affected in the turn and direction.

4 The Lens Stop & Table:

The spectacles & holder where the spectacles are held where a holder is replaced the wheel for decentric and centric alignment.

3 The Power Drum:

The lower drum is a graduated drum. All labels within can only that the power marking would carry as shown and its accurate per cent.




Decoding the Reticle Target


Step 1: Unneutralized → **Step 2: Target 1 (Sphere)** → **Step 3: Target 2 (Cylinder)**



Target is a blurry array of lines.



The three thin lines come into sharp, glowing cyan focus. The power drum reading at this exact moment dictates the sphere power.




The three thick lines come into sharp focus. The dioptric difference between Target 1 and Target 2 dictates the cylinder power.



CLINICAL RULE: Always neutralize Target 1 (thin lines) before Target 2 (thick lines) when working in minus cylinder format.

Calculating the Presbyopic Add


Translating (Lined Bifocal)



Process: Measure the primary distance power first. Move the lens table *up* to position the near segment over the lens stop. Read the new power.

Formula: The dioptric difference between distance and near = The Add.

Aspheric (Progressive Addition Lens)



Process: Locate the temporal circular engraving to find the pre-stamped Add power, or utilize the automated lensometer's dedicated PAL mode.

Identifying Prism Displacement


Each concentric ring represents exactly 1 diopter of prism.

Prism is simply the optical center moving away from the physical center.



Diagnostic: In this specific visual example, the target displacement toward the top right indicates Base Up and Base Out prism.

Determining Refractive Index with the Lens Clock



Step 1: Read total power on lensometer (e.g., $-5.00D$).

Step 2: Measure Front Curve (+) and Back Curve (-) with the clock. Add them together.

Step 3: Cycle the clock's index setting (1.49, 1.53, 1.56, 1.59, 1.61, 1.67, 1.71, 1.74) until the calculated curve matches the lensometer power.

Example: At index 1.61, Front reads $+0.50D$ and Back reads $-5.50D$. Total = $-5.00D$. The material index is 1.61.

The Technician's Lensometry Workflow

